


DUAL-CURED CORE BUILD-UP MATERIAL

CLEARFIL™ DC CORE PLUS




SIMPLE! STRONG! CONFIDENT! The Core Build-Up Material.

CLEARFIL™ DC CORE PLUS is a dual-cure, radiopaque two-component core build-up material supplied in a convenient automix delivery system.

CLEARFIL™ DC CORE PLUS provides the dentist with the confidence of durable, lasting restorations.

CLEARFIL™ DC CORE PLUS has superior compressive and flexural strength with a filler loading of 74wt%(52vol%). **CLEARFIL™ DC CORE PLUS** cuts like dentin when being prepared, and affords long-term durability.

A strong core delivers the confidence of a great restoration.



Product Range: CLEARFIL™ DC CORE PLUS


ORDERING INFORMATION:

CLEARFIL™ DC CORE PLUS

#2942-WD Dentin
#2943-WD White

Content:

- Automix Syringe (1pc)
 - A Paste (8.9 g, 4.5 ml)
 - B Paste (9.0 g, 4.5 ml)
- Mixing Tips (20pcs)
- Guide Tips (L) (10pcs)
- Guide Tips (S) (10pcs)



Dentin



White

ACCESSORIES

#1029-WD

CLEARFIL™ DC CORE PLUS
Guide Tip (L) 20 pcs



#1031-WD

CLEARFIL™ DC CORE PLUS
Mixing Tip 20 pcs



#1030-WD

CLEARFIL™ DC CORE PLUS
Guide Tip (S) 20 pcs



- Read instructions for use supplied with the product before use.
- The measurement data contained in this brochure are obtained by the results of tests performed by Kuraray Noritake Dental Inc.
- The colors of the product pastes and hardened products shown in this brochure may vary from the actual colors.
- The specifications and appearance of the product are subject to change without notice.

Kuraray Noritake Dental Inc.

Ote Center Bldg.
1-1-3, Otemachi, Chiyoda-Ku, Tokyo 100-0004, Japan
Website : <http://www.kuraraynoritake.com/wd>